





HAZARD MANAGEMENT – SAFE OPERATING PROCEDURE (SOP)

Only to be completed where required as a control measure under a Risk Assessment

NAME OF THE TASK/ACTIVITY	OLYMPUS FV3000 CONFOCAL MICROSCOPE	DATE: 24/02/2020
LOCATION	ADELAIDE MICROSCOPY, MLS BUILDING RM 1.22 AND AHMS LEVEL 6	Insert photo (Optional)
RISK ASSESSMENT (RA) NAME	Olympus FV3000 Confocal Microscope	
Residual risk rating on the RA	<input checked="" type="checkbox"/> Low <input type="checkbox"/> Medium <input type="checkbox"/> High <input type="checkbox"/> Very High	
Hazards identified on the RA	Contact with electricity or potential for electric shock Exposure of eyes to the Laser beam. Exposure to non-ionising radiation (UV light) Mercury lamp breakage and vapour or liquid leak	
PERSONAL PROTECTIVE EQUIPMENT		
	Eye protection: <input checked="" type="checkbox"/> Safety glasses <input type="checkbox"/> Eye shields <input type="checkbox"/> Safety goggles <input type="checkbox"/> Other:	
	Hand protection: <input type="checkbox"/> Rubber <input type="checkbox"/> Cut resistant <input type="checkbox"/> Leather <input type="checkbox"/> Vinyl <input type="checkbox"/> Neoprene <input checked="" type="checkbox"/> Nitrile <input type="checkbox"/> Barrier creams <input type="checkbox"/> Other:	
	<input checked="" type="checkbox"/> Enclosed footwear: <input type="checkbox"/> Footwear that is resistant to spills of hazardous substances <input type="checkbox"/> Boots with steel caps <input type="checkbox"/> Other:	
	Protective clothing: <input checked="" type="checkbox"/> Lab coat <input type="checkbox"/> Gown <input type="checkbox"/> Long sleeves <input type="checkbox"/> Long pants <input type="checkbox"/> High visibility <input type="checkbox"/> Helmet <input type="checkbox"/> Sun protection <input type="checkbox"/> Other:	
DESCRIBE, IN SEQUENCE, STEPS TO COMPLETE THE ACTIVITY SAFELY		
Pre-operational checks		
YOU MUST NOT USE THIS MACHINE UNTIL YOU HAVE HAD APPROPRIATE TRAINING BY TRAINED ADELAIDE MICROSCOPY STAFF. Unauthorised use may result in damage to the instrument.		
Operational checks/steps to complete the activity from start to finish (including transport and waste disposal where relevant)		
General		
The Olympus FV3000 scanning confocal is a fluorescence microscope in which lasers are used to image thin optical sections of biological samples. Dyes in the samples fluoresce and the emitted light passes through a very fine pinhole to detectors. Detectors are connected to a computer which builds up the image. UV light is also used to locate specimens using wide-field fluorescence.		
The Olympus FV3000 confocal uses both UV light and lasers as sources of illumination. The instrument utilises Class3B Lasers, including solid state 405, 445, 488, 514, 561, 594 and 640 lasers. UV light is provided by an Olympus U-HGLGPS.		
Hazards		
Potential for electric shock if a user were to remove panels from the microscope.		
Exposure of eyes to the Laser beam can cause blindness.		
Exposure of eyes to UV light can cause eye damage.		
Breakage of the hot or cold lamp may release mercury vapour or liquid.		
Risk Control Measures		
<u>Engineering controls:</u>		
The user operable parts on the Olympus FV3000 scanning confocal microscope are all accessible from the front of the instrument, and include the inverted light microscope and controller, the automated stage controller, and the computer (switch, mouse and keyboard). There is no risk involved in the operation of these parts. However, misuse of these parts can result in damage to the instrument. Users of the instrument should not remove any fixture or panel from the microscope or access the rear of the instrument.		
The microscope has several inbuilt safety features to stop people injuring themselves with the laser or the mercury lamp.		

HSW Handbook	Hazard Management	Effective Date:	17 December 2019	Version 3.0
Authorised by	Chief Operating Officer (University Operations)	Review Date:	17 December 2022	Page 1 of 2
Warning	This process is uncontrolled when printed. The current version of this document is available on the HSW Website.			

The laser beams are contained and shielded in boxes which cannot be opened accidentally. An Interlock must be manually switched off in the computer software before the lasers can be used; the Interlock prevents laser light entering the binocular eyepieces. When the lasers are in operation it is only possible to view the sample on the computer screen (it is not possible to view the sample down the microscope binocular).

There are low risks to the user associated with use of the U-HGLGPS UV lamp. The microscope has several inbuilt safety features to stop people injuring themselves with the UV lamp.

The Olympus U-HGLGPS lamp is pre-centered and requires no alignment, reducing the risk of overheating and explosion.

The UV light source is contained and shielded in a box which cannot be opened accidentally. To allow the UV light to the specimen, the operator must manually operate a button on the computer software and/or a shutter on the Intensilight controller. The UV light may also be turned on and off via a shutter on the microscope.

An external orange filter guards the user's eyes from the stage area. The only way to be exposed to UV or laser light at the microscope stage is to deliberately remove the guard or reach around behind the guard.

The lifetime of the Olympus U-HGLGPS pre-centered mercury lamp is 2,000 hours, which is monitored by a counter on the front of the lamp housing. There is no risk of incorrectly aligning a pre-centred lamp. The lamp does not have overpressure when cool.

Procedural controls:

Only trained users to operate the instrument. All new users are to be given practical training in instrument operation by a member or Adelaide Microscopy staff. Users must also follow guidelines in the manual and safe operating procedures for operation of the microscope and UV light source.

Users must not tamper with the instrument such that they can look directly into a laser beam or a reflection of a laser beam. Users must not introduce any reflective objects into the laser beam path – e.g. micromanipulators

Handling of biological material may present hazards; the safe operating procedures for handling biological material must be followed. The handling of other laboratory items (for example, sharps, clearing agents and chemicals) must follow the relevant safe operating procedures.

General Procedures:

Users should operate the instrument in accordance with the manufacturer supplied operating instructions under the instruction of a member of Adelaide Microscopy staff.


On completion of work – steps to make safe (including clean up, any waste disposal & service/maintenance requirements)

Follow the shutdown procedure in the manual.

Emergency and Spill Procedures, Transport or storage requirements (where relevant), First aid/Medical

In the event of an injury, please advise an Adelaide Microscopy staff member and first aid officer for treatment and the local HSW representative to report the incident.

Prepared by

People involved in the drafting of this SOP	Jane Sibbons Aoife McFadden		
Person authorising the SOP	Name:	Angus Netting	Signature 
	Position:	Director, Adelaide Microscopy	
This SOP must be reviewed after any incident/injury associated with this activity or when a Risk assessment is reviewed.			

HSW Handbook	Hazard Management	Effective Date:	17 December 2019	Version 3.0
Authorised by	Chief Operating Officer (University Operations)	Review Date:	17 December 2022	Page 2 of 2
Warning	This process is uncontrolled when printed. The current version of this document is available on the HSW Website.			