INDEX, ÍNDICE, EYPETHPIO

1-11	ENGLISH
12-23	ESPAÑOL
24-35	ITALIANO
36-47	DEUTSCHE
48-59	ΕΛΛΗΝΙΚΟ
60-71	FRANÇAIS
72-73	SYMBOLS SIMBOLOS SIMBOLI SYMBOL SYMBOLES ΣΥΜΒΟΛΟ SIMBOLIS SYMBOLY

SENSITITRE®
YEASTONE®
For *in vitro* Diagnostic Use

For full plate information, including plate layout, QC information please refer to www.trekds.com/techinfo. The plate code and batch number will be required.

INTENDED USE

The Sensititre® susceptibility system is an *in vitro* diagnostic product for susceptibility testing of non-fastidious yeast including *Candida* species, *Cryptococcus* species, *Aspergillus* species and miscellaneous other rapid growing yeast species.

It is a micro broth method that provides qualitative quantitative Minimum Inhibitory Concentration (MIC) results in a dried plate format. TREK Diagnostic Systems manufactured broth has only been validated with Sensititre® Products.

PRINCIPLES OF USE

The Sensititre® yeast susceptibility test is a colorimetric microdilution test. Each plate is dosed with antifungal agents at appropriate dilutions, and a colorimetric indicator.

Results are read manually by observing the lowest antifungal concentration showing inhibition of growth (as evidenced by no color change).

PRECAUTIONS

Results should be used as an aid to selecting the drug of choice for treatment

Only personnel trained in susceptibility testing techniques should use the system

Since living microorganisms used with this product can be infectious to the user, proper handling and disposal methods should be used.

STORAGE AND SHELF LIFE

The plates should be stored at room temperature (15-25°C) away from direct sunlight and direct heat. Each plate is packaged in foil with a silica gel desiccant. Do not use the plate or broth if past its expiration date, or the desiccant color is not blue or orange or the foil pouch is damaged.

Inoculate the plate within 5 hours of removal from the pouch.

PROCEDURE

Materials Included:

YeastOne® susceptibility plate

Adhesive seal

Materials Not Included [TREK Inc Product code]:

Sensititre® demineralised water [T3339]

Sensititre® yeast susceptibility inoculum broth [Y3462]

Sensititre® doseheads (for use with AutoInoculator®/AIM®) [E3010]

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Sensititre® Vizion® [V2021]

Sensititre® Nephelometer® [V3011]

Manual Viewer [V4007]
0.5 McFarland turbidity standard [E1041]
Bacteriological loop
20 μ l pipette
Sterile inoculum reservoir
100 μ l pipettor and disposable tips
Quality control strains
Fungal growth medium agar plates e.g Sabauroid dextrose agar (SDA)
Incubator 34-36°C, non-CO₂
Vortex mixer
Current CLSI documents
Aspergillus spp only
Takashio or potato dextrose agar
Cotton Swab
Sterile saline with 0.05% Tween 20
Spectrophotometer

SPECIMEN COLLECTION AND PREPARATION

Specimens should be collected, transported, stored and then plated on to primary isolation medium to give isolated colonies using standard procedures. (1)

SELECTION OF SUSCEPTIBILITY TEST BROTH

Sensititre® approved broths are performance tested for use in Sensititre® susceptibility products.

INOCULATION PROCEDURE (*Candida* and *Cryptococcus* spp.)

Allow all broths to come up to room temperature before use.

A final organism density of approximately 1.5 - 8 X 10³ CFU/ml is recommended.

1. Pick several well-isolated colonies of >1mm diameter from a pure 24-hour culture of the Yeast isolate, and emulsify into sterile water. Vortex mix the suspension for 15 seconds, ensuring that the suspension is uniform. If clumping occurs, allow them to settle before adjusting the density. Adjust to a 0.5 McFarland standard visually or using a Sensititre® Nephelometer®
 2. Transfer 20 μ l of the suspension into 11 ml of YeastOne® inoculum broth, to give an inoculum of 1.5 - 8 X 10³ CFU/ml.
- Steps 1 and 2 should be completed in 15 minutes.
3. Transfer 100 μ l by either:
 - a. **Sensititre® AutoInoculator®/ Sensititre AIM®.** Replace the tube cap with a Sensititre® single-use dosehead and inoculate the plate according to the AutoInoculator® /AIM® instructions.

Remove the test tube dosehead combination from the AutoInoculator® /AIM® within 30 seconds of dosing a plate and store inverted in a rack or discard.

- b. **Manual pipette.** Pour the broth into a sterile seed trough and inoculate the plate using an appropriate pipette.

Inoculate broth into a plate within 15 minutes.

4. A check of the colony count should be done by removing 10 μ l from the positive control well and plating onto Sabauroid dextrose agar (SDA). A correct inoculum will produce 10-80 colonies.

5. Cover all wells with the adhesive seal. Avoid creases as these can lead to skips.

INOCULATION PROCEDURE (*Aspergillus* spp.) (5-7, 22)

Allow all broths to come up to room temperature before use.

1. Subculture from Sabouraud dextrose agar onto Tasashio or potato dextrose agar. Incubate for 7 days at 35°C to obtain adequate sporulation.
 2. Collect conidia with a cotton swab and suspend in sterile saline with Tween.
 3. Allow heavy particles to settle for 3 to 5 minutes.
 4. Collect supernatant and mix with a vortex mixer.
 5. Adjust turbidity of supernatant to 80 to 82% transmittance at 530nm measured with a spectrophotometer equivalent to an inoculum of $0.6 - 5 \times 10^6$ cfu/ml. Alternatively, adjust to a 0.5 McFarland Standard (22).
 6. Add 100µl to 11ml of YeastOne® inoculum broth to give a final inoculum of $0.5-5 \times 10^4$ cfu/ml.
 7. Transfer 100µl by either:
 - a. **Sensititre® AutoInoculator®/Sensititre® AIM®**. Replace the tube cap with a Sensititre® single-use dosehead and inoculate the plate according to the AutoInoculator®/AIM® instructions.
- Remove the test tube dosehead combination from the AutoInoculator® /AIM® within 30 seconds of dosing a plate and store inverted in a rack or discard.
- b. **Manual pipette**. Pour the broth into a sterile seed trough and inoculate the plate using an appropriate pipette.
Inoculate broth into a plate within 15 minutes.
 8. A check of the colony count should be done by removing 10µl from the positive control well and plating onto Sabauroid dextrose agar (SDA). A correct inoculum will produce 50-500 colonies.
 9. Cover all wells with the adhesive seal. Avoid creases as these can lead to skips.

INCUBATION

Minimally incubate the plates for 24-25 hours at 35°C in a non-CO₂ incubator.

Cryptococcus species should be incubated for 72 hours.

Aspergillus species should be incubated for 48 to 72 hours.

Incubation at temperatures over 35°C may affect the performance of these plates.

Up to 3 plates can be stacked if not incubated in the ARIS®.

READING TEST RESULTS

Plates may be read visually under normal laboratory lighting using a manual mirror viewer or by using the Sensititre Vizion® System. Refer to the Vizion User manual for additional information. Yeast growth in the antifungal solutions will be evident as a change in the colorimetric growth indicator from blue (negative) to red (positive). Some yeasts may not change the indicator completely to red, but display more of a purpling of the indicator. Some organisms may show a slight purpling in posaconazole, voriconazole, fluconazole, itraconazole and ketoconazole. (see details for reading below).

1. Examine the positive growth well after 24 hours incubation (*Candida* species). If the growth well is red, the endpoints for the antifungals can be interpreted. If the well is still blue or only faintly purple, re-incubate for an additional 24 hours and re-examine.

DO NOT READ TURBIDITY IN THE SENSITITRE YEASTONE® PLATES.**Read Only Color Change.**

2. The MIC is the lowest concentration of an antifungal agent that substantially inhibits growth of the organism as detected by a color change. The amount of color change in the wells containing the agent should be compared to the color of the positive growth-control wells
3. No growth occurs when there is no change in the blue indicator in any dilution of an antifungal. The organism is susceptible to the lowest concentration of antifungal
4. The MIC is recorded as the lowest concentration of antifungal agent preventing the development of a red or purple growth well, i.e. first blue. .
5. The organism is resistant to the highest concentration of antifungal when growth is seen in all wells. The MIC endpoint should be recorded as "greater than" (>) the highest concentration.
6. For *Aspergillus* species read the MIC as the lowest concentration with a blue color.

INTERPRETATION OF RESULTS**TABLE 2. Illustration and the interpretation of test results that may occur**

Well Concentration µg/ml	1	2	4	8	16	32	R = RED: Positive growth indication B = BLUE: Negative growth indication
A.	R	R	R	B	B	B	Typical growth pattern; MIC endpoint is 8 µg/ml.
B.	R	R	R	R	R	R	Growth in all wells; MIC endpoint is >32 µg/ml.
C.	B	B	B	B	B	B	No growth in any well; MIC endpoint is ≤1 µg/ml.
D.	R	R	R	B	R	R	"Skipped Well". MIC endpoint is >32 µg/ml. Disregard "skip" when wells on either side have growth. If more than one "skip" should occur in a column, the test results are invalidated ¹
E.	R	R	B	B	R	R	Double "Skipped Well". The test should be repeated ¹

¹ With careful technique these occurrences are uncommon.

READING NOTES


Amphotericin B. For amphotericin B at 24 hours, the endpoints are typically easily defined and the MIC is read as the lowest drug concentration that prevents any discernible color change. Trailing endpoints with Amphotericin B are not usually encountered.

The first well showing a distinct color change as compared to the positive growth well is the MIC.

Flucytosine and Azole Antifungals. *Candida albicans*, *C. glabrata* and *C. tropicalis* with flucytosine and azoles, such as fluconazole, itraconazole, ketoconazole, voriconazole and posaconazole may give endpoints that are typically less sharp because of trailing growth, and may be a significant source of variability. Trailing occurs when a slight colour change persists and it is often identical for all drug concentrations above the MIC. The MIC should be read as the first well showing a less intense colour change compared to the positive growth control well. Reference strains of defined susceptibility may also help to train personnel. Isolates of *Candida krusei* are assumed to be intrinsically resistant to fluconazole and their MICS should not be interpreted, (1) a comment should accompany the test result reported.

Trailing endpoint: This occurs when a slight color change persists and is often identical in several concentrations. The MIC should be read as the first well showing a less intense color change compared to the positive growth control well.

Echinocandins. The MIC end points should be determined after 24 hours of incubation at 35°C. The MIC should be read as the first well showing a less intense color change as compared to the positive control well.

Itraconazole:


Itraconazole can occasionally come out of solution at concentrations of $\geq 4 \mu\text{g/ml}$. This can result in the affected well exhibiting growth and turning red.

From time to time we are encountering paradoxical growth in the higher concentrations of Itraconazole on the Sensititre yeast panels which results in pinking of these wells.


The paradoxical effect also known as the Eagle phenomenon refers to an observation where an increase in the antimicrobial concentration beyond a certain point paradoxically

results in an increase to the number of bacteria that survive. An explanation could be that as the concentration is too high, the agent might be self-antagonising the receptor with which it binds (penicillin binding proteins, for example, in the case of a penicillin).

Resolution

The growth in the high concentration should be ignored unless you have growth in all of the other concentrations of Itraconazole. In the example below the 0.5μg/ml well highlighted with the black square is where the MIC result should be recorded.

If you have any other questions or concerns please contact the technical support department on +441342318777 or email techsupport@trekds.co.uk.

Alternatively please contact +1-800-871-8909 or email techsupport@trekds.com in the USA.

Contamination/ Skips

Alternatively, a pink (growth) well between blue (no growth) wells could be indicative of contamination. Sub-culture well contents to ascertain the cause.

A blue well in a series of red growth wells indicates a “skip” and should be ignored. The MIC should be read above any skip wells. If there is more than one skipped well, the antifungal should not be reported.

QUALITY CONTROL

Frequency of quality control testing should be according to local guidelines (1). Inoculum should be cultured onto a suitable medium to check for purity. Test results are invalid if a mixed culture is detected.

All Sensititre® plates include positive control wells. Tests are invalid unless there is distinct growth in all positive control wells

For user quality control of the MIC system, the following culture(s) from the American Type Culture Collection (ATCC®) are recommended:

<i>Candida krusei</i> *	ATCC® 6258
<i>Candida parapsilosis</i>	ATCC® 22019

*ATCC now lists these organisms as *Issatchenkia orientalis*.

Results should **not** be reported if QC results are not in range.

Expected QC values are provided in Table 3

Contact Sensititre® Distributor or TREK Diagnostic Systems in the event that quality control discrepancies cannot be resolved.

TABLE 3. Recommended 24 and 48-hour MIC ranges ($\mu\text{g/ml}$) or Quality Control strains.

Ranges that are different from CLSI quality control ranges are underlined (1).

Antifungal Agent	<i>Candida krusei</i> ATCC 6258		<i>Candida parapsilosis</i> ATCC 22019	
	24 hour	48 hour	24 hour	48 hour
Amphotericin B	0.5 - 2	1 - 4	0.25 - 2	0.5 - 4
Anidulafungin	0.03 - 0.12	-	0.25 - 2	-
Caspofungin	0.12 - 1	0.25 - 1	0.25 - 1	0.5 - 4
Fluconazole	8 - 64	16 - 128	0.5-4	<u>2 - 8</u>
5-Flucytosine	4 - 16	8 - 32	<u>0.12 - 0.5</u>	0.12 - 0.5
Itraconazole	0.12 - 1	0.25 - 1	0.12 - 0.5	0.12 - 0.5
Ketoconazole	0.12 - 1	0.25 - 1	0.03 - 0.25	0.06 - 0.5
Micafungin	<u>0.06 - 0.25</u>	0.12 - 0.5	0.5 - 2	0.5 - 4
Posaconazole	0.06 - 0.5	0.12-1	0.06 - 0.25	0.06 - 0.25
Voriconazole	0.06 - 0.5	0.12 - 1	0.015 - 0.12	0.03 - 0.25

TABLE 4. MIC Interpretative Criteria ($\mu\text{g/ml}$) for *Candida* Species as per CLSI M27

Antifungal Agent	Susceptible	Susceptible dose dependent	Intermediate	Resistant	Non-susceptible
Anidulafungin	<u>≤ 2</u>				> 2
Caspofungin	<u>≤ 2</u>				> 2
Fluconazole*	<u>≤ 8</u>	16 - 32		≥ 64	
5-Flucytosine	<u>≤ 4</u>		8 - 16	≥ 32	
Itraconazole	<u>≤ 0.125</u>	0.25 - 0.5		≥ 1	
Micafungin	<u>≤ 2</u>				> 2
Voriconazole	<u>≤ 1</u>	2		≥ 4	

Isolates of *Candida krusei* are assumed to be intrinsically resistant to fluconazole and their MICs should not be interpreted using this scale

NOTE 1: Shown are the breakpoints ($\mu\text{g/mL}$) for *Candida* spp. Against the indicated agents. If minimal inhibitory concentrations (MICs) are measured using a scale that yields results falling between categories, the next higher category is implied. Thus an isolate with a fluconazole MIC of 12.5 $\mu\text{g/mL}$ would be placed in the S-DD category.

Please refer to CLSI (1) for more information concerning interpretation of results.

LIMITATIONS

1. Sensititre YeastOne® plates are for use with non-fastidious yeast including *Candida* species, *Cryptococcus* species, *Aspergillus* species and miscellaneous rapid growing

yeast species. They are not intended for fastidious or slow growing yeast such as *Histoplasma* or *Blastomyces*, and filamentous fungi.

2. Comparison between the Sensititre YeastOne® at 24 hours and the CLSI reference method at 48 hour was evaluated. However due to the difficulty in correlating end points of trailing organisms (*C. albicans*) at 48 hours incubation, high error rates are observed.

3. Testing of fungi and antifungal agents is inherently less precise than testing bacteria.

4. Some investigators believe the 24-hour reading is more appropriate than the 48-hour reading because of the problem with trailing with certain isolates. The CLSI official standard indicates that readings should be accomplished at 48 hours. Until sufficient data is collected and analyzed, the question of most clinically relevant time of reading remains unanswered. Reporting of results should indicate clearly the times of reading.

5. For additional guidance, review of CLSI Antifungal Susceptibilities Standard M27 is encouraged.

6. Colour change is the indicator of the end point, not turbidity. (This fact alleviates some major concerns with the interpretation of certain *Candida* species because of 'trailing'. Trailing is more commonly seen with isolates other than those of blood and other sterile body fluids.)

7. Do not read at 24 hours if the control well has not completely turned positive.

8. The performance of voriconazole with *Cryptococcus* species and rapid growing yeast species has not been determined. Voriconazole MIC's should therefore only be reported for *Candida* species.

9. Use only with Sensititre® approved yeast susceptibility inoculum broth. The use of other broths could result in error.

10. As with any in-vitro susceptibility testing method, the results of testing should be correlated with the patient's clinical response to prescribed therapy.

11. Performance has only been established for Amphotericin B, Itraconazole Posaconazole and Voriconazole with *Aspergillus* species.

12. Reference 5 showed high level (>99% agreement) with CLSI method for amphotericin B and *Aspergillus* species. Lower agreement was seen with itraconazole, *A. fumigatus*, *A. flavus* and *A. terreus* had >90% agreement whilst *A. nidulans* had 85% and *A. ustus* 33% on panels incubated for 48 hours with a 10^3 cfu/ml inoculum. Agreement over 90% was observed with an inoculum of 10^4 cfu/ml and 72 hours incubation with itraconazole and amphotericin B (6).
13. Correlation of the MIC for Caspofungin to the treatment outcome following caspofungin use has not been fully established. (9)
14. Performance of the YeastOne® with Anidulafungin, Caspofungin, and Micafungin with Cryptococcus, Aspergillus species and rapid growing yeast species **other than Candida species** has not been established. Anidulafungin, Caspofungin and Micafungin MIC's should therefore only be reported for Candida species.
15. Performance of the YeastOne® with Posaconazole with Cryptococcus has not been established.
16. Only instruments supported by Sensititre i.e. a simple mirror viewer, Sensitouch, Vizion, Sensititre Autoreader , Optiread and ARIS instruments must be used to report results with CE IVD and FDA cleared Sensititre products, any other system used will not be supported.

PERFORMANCE

Panels read either manually or are designed to give comparable performance to CLSI reference microbroth procedure. Comparable performance is defined as \geq 90% agreement to within a doubling dilution of the reference MIC (1).

For further information contact TREK Diagnostic systems or your local distributor

BIBLIOGRAPHY

- (1) Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087.
- (2) Espinel-Ingroff, A., Pfaller, M., Messer, S.A., Knapp, C.C., Killian, S., Norris, H.A., Ghannoun M.A. (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS, M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology*. **37**: 591-595.
- (3) Espinel-Ingroff, A., Knapp, C.C., Holiday, N., Killian, S. (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection*. **8**: Suppl 1. Abstract 0125.
- (4) Davey, K.G., Szekely, A., Johnson, E.M., Warnock,D.W., (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Meletiadis, J., Moutin, J.W., Meis, F.G.M., Bouman, B.A., Verweij, P.E., and EuroFung Network. (2002). Comparison of the E Test and the Sensititre® Colorimetric Methods with the NCCLS, Proposed Standard for Antifungal Susceptibility Testing of *Aspergillus* species. *Journal of Clinical Microbiology*. **40**: 2876-2885.

- (6) Espinel-Ingroff, A., et al. (1997). Multicenter Evaluation of Proposed Standardized Procedure for Antifungal Susceptibility Testing of Filamentous Fungi. *Journal of Clinical Microbiology*. **35**: 139-143.
- (7) Pemán, J., et al. (2001). Comparison of the Sensititre YeastOne® Colorimetric method against reference M38-P for testing susceptibility of *Aspergillus* to Amphotericin B and Itraconazole. *Clinical Microbiology and Infection*, **7**: Suppl. 1. Abstract P694.
- (8) Linares, M., J. G. Charriel, F. Solís, F. Rodríguez, A. Ibarra and M. Casal. (2005) Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253
- (9) Kartsonis, N., et al. (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623
- (10) Pfaller, M. A., et al. (2004) Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases* **50** :113-117
- (11) Pfaller, M. A., A. Espinel-Ingroff and R.N. Jones. (2004) Clinical evaluation of the Sensititre YeastOne® antifungal susceptibility testing of the new triazoles voriconazole, posaconazole and ravuconazole. *Journal of Clinical Microbiology* **42**: 4577-4580
- (12) Espinel-Ingroff, A., et al. (2004) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the NCCLS, M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology* **42**: 718- 721
- (13) Linares, M., J., G. Charriel, F. Solís and M. Casal. (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology* **42**: 899-902
- (14) Espinel-Ingroff, A., M. Pfaller, S.A. Messer, C.C. Knapp, S. Killian, H.A. Norris, and M.A. Ghannoum. (1999) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology* **37**: 591-595
- (15) Barry, A.L., M.A. Pfaller, S.D. Brown, A. Espinel-Ingroff, M.A. Ghannoum, C.c. Knapp, R.P. Rennie, J.H. Rex, and M.G. Rinaldi. (2000) Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology* **38**: 3457-3459
- (16) Canton, E., et al (2005) .Sensititre YeastOne® Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of NCCLS, M27-A2 multicenter study . *Antimicrobial Agents and Chemotherapy*. **49**: 1604 -1607
- (17) Holiday, N.M..., et al (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne® susceptibility plate. *American Microbiology Abstract C-191*

- (18) Torres-Nabona, M., J. Guinea, J.Martinez_Alarcon, T. Pelaez, and E. Bouza. (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of NCCLS, M38-A , Sensititre_YeastOne® and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 - 1129
- (19) Espinel-Infogroff, A. (2006) Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622
- (20) Patel, R., C. Mendrick, C.C. Knapp, R. Grist, and P.M. McNicholas. (2007) Clinical evaluation of the Sensititre YeastOne® plate for testing susceptibility of filamentous fungi to posaconazole. *Journal of Clinical Microbiology* **45**: 2000-2001
- (21) Alexander, B. D., T.C. Byrne, K.L. Smith, K.E. Hansen, K.J. Anstrom, J.R. Perfect and L. B. Reller. (2007) Comparative evaluation of Etest and Sensititre YeastOne® panels against the Clinical Laboratory Standards Institute M27-A2 reference broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology* **45**: 698-706
- (22). Patel, R., *et al.* (2007). Clinical evaluation of the Sensititre YeastOne plate for testing susceptibility of filamentous fungi to posiconazole. *Journal of Clinical Microbiology*. **45**: 2000-2001

DISCLAIMER

The information provided in this technical insert is current at the time of printing and may change without notice.

The latest information can be downloaded from www.trekds.com\techinfo or by Contacting TREK Technical services.

Manufactured by TREK Diagnostic Systems
Units 17-19 Birches Industrial Estate, East Grinstead,
West Sussex. RH19 1XZ, UK.
Tel: +44- 1342-318777

Distributed by TREK Diagnostic Systems,
982 Keynote Circle, Suite 6, Cleveland, Ohio 44131
Technical Service USA: 1 (800) 642-7029

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- (1) Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. M27 Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087.
- (2) Espinel-Ingroff, A., Pfaller, M., Messer, S.A., Knapp, C.C., Killian, S., Norris, H.A., Ghannoun M.A. (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the National Committee for Clinical Laboratory Standards M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology* **37**: 591-595.
- (3) Espinel-Ingroff, A., Knapp, C.C., Holiday, N., Killian, S. (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection* **8**: Suppl 1. Abstract 0125.
- (4) Davey, K.G., Szekely, A., Johnson, E.M., Warnock,D.W., (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Meletiadis, J., Moutin, J.W., Meis, F.G.M., Bouman, B.A., Verweij, P.E., and EuroFung Network. (2002). Comparison of the E Test and the Sensititre® Colorimetric Methods with the NCCLS Proposed Standard for Antifungal Susceptibility Testing of *Aspergillus* species. *Journal of Clinical Microbiology*. **40**: 2876-2885.
- (6) Espinel-Ingroff, A., etal. (1997). Multicenter Evaluation of Proposed Standardized Procedure for Antifungal Susceptibility Testing of Filamentous Fungi. *Journal of Clinical Microbiology* **35**: 139-143.
- (7) Pemán, J., etal. (2001). Comparison of the Sensititre YeastOne® Colorimetric method against reference M38-P for testing susceptibility of *Aspergillus* to Amphotericin B and Itraconazole. *Clinical Microbiology and Infection*, **7**: Suppl. 1. Abstract P694.

- (8) Linares,M, J., G. Charriel, F. Solís, F. Rodriguez, A. Ibarra and M. Casal. (2005) Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253
- (9) Kartsonis, N., etal. (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623
- (10) Pfaller, M, A., etal. (2004) Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases* **50** :113-117
- (11). Pfaller, M, A., A. Espinel-Ingroff and R.N. Jones. (2004) Clinical evaluation of the Sensititre YeastOne® antifungal susceptibility testing of the new triazoles voriconazole, posaconazole and ravuconazole. *Journal of Clinical Microbiology* **42**: 4577-4580
- (12). Espinel-Ingroff, A., etal. (2004) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the NCCLS M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology* **42**: 718- 721
- (13) Linares,M, J., G. Charriel, F. Solís and M. Casal. (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology* **42**: 899-902
- (14) Espinel-Ingroff, A., M. Pfaller, S.A. Messer, C.C. Knapp, S. Killian, H.A. Norris, and M.A. Ghannoum. (1999) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology* **37**: 591-595
- (15) Barry, A.L., M.A. Pfaller, S.D. Brown, A. Espinel-Ingroff, M.A. Ghannoum, C.c. Knapp, R.P. Rennie, J.H. Rex, and M.G. Rinaldi. (2000) Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology* **38**: 3457-3459
- (16) Canton, E., etal (2005) .Sensititre YeastOne® Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of NCCLS M27-A2 multicenter study . *Antimicrobial Agents and Chemotherapy*. **49**: 1604 -1607
- (17) Holiday, N.M.., etal (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne® susceptibility plate. *American Microbiology Abstract C-191*
- (18) Torres-Nabona, M., J. Guinea, J.Martinez_Alarcon, T. Pelaez, and E. Bouza. (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of CLSI M38-A , Sensititre_YeastOne® and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 - 112

- (19) Espinel-Infgriff, A. (2006) Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622
- (20) Patel, R., C. Mendrick, C.C. Knapp, R. Grist, and P.M. McNicholas. (2007) Clinical evaluation of the Sensititre YeastOne® plate for testing susceptibility of filamentous fungi to posaconazole. *Journal of Clinical Microbiology* **45**: 2000-2001
- (21) Alexander, B. D., T.C. Byrne, K.L. Smith, K.E. Hansen, K.J. Anstrom, J.R. Perfect and L. B. Reller. (2007) Comparative evaluation of Etest and Sensititre YeastOne® panels against the Clinical Laboratory Standards Institute M27-A2 reference broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology* **45**: 698-706
- (22). Patel, R., *etal.* (2007). Clinical evaluation of the Sensititre YeastOne plate for testing susceptibility of filamentous fungi to posiconazole. *Journal of Clinical Microbiology*. **45**: 2000-2001

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029-Yeast-ROW-IVD-E-V2.0
Date Updated: 06 August 2012

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BIBLIOGRAFIA

- (1) Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. M27 Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087.
- (2) Espinel-Ingroff, A., Pfaller, M., Messer, S.A., Knapp, C.C., Killian, S., Norris, H.A., Ghannoun M.A. (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology* **37**: 591-595.
- (3) Espinel-Ingroff, A., Knapp, C.C., Holiday, N., Killian, S. (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection* **8**: Suppl 1. Abstract 0125.
- (4) Davey, K.G., Szekely, A., Johnson, E.M., Warnock,D.W., (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Meletiadis, J., Moutin, J.W., Meis, F.G.M., Bouman, B.A., Verweij, P.E., and EuroFung Network. (2002). Comparison of the E Test and the Sensititre® Colorimetric Methods with the NCCLS Proposed Standard for Antifungal Susceptibility Testing of *Aspergillus* species. *Journal of Clinical Microbiology*. **40**: 2876-2885.
- (6) Espinel-Ingroff, A., etal. (1997). Multicenter Evaluation of Proposed Standardized Procedure for Antifungal Susceptibility Testing of Filamentous Fungi. *Journal of Clinical Microbiology* **35**: 139-143.
- (7) Pemán, J., etal. (2001). Comparison of the Sensititre YeastOne® Colorimetric method against reference M38-P for testing susceptibility of *Aspergillus* to Amphotericin B and Itraconazole. *Clinical Microbiology and Infection*, **7**: Suppl. 1. Abstract P694.
- (8) Linares,M, J., G. Charriel, F. Solís, F. Rodriguez, A. Ibarra and M. Casal. (2005) Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253
- (9) Kartsonis, N., etal. (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623

- (10) Pfaller, M. A., *et al.* (2004) Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases* **50** :113-117
- (11). Pfaller, M. A., A. Espinel-Ingroff and R.N. Jones. (2004) Clinical evaluation of the Sensititre YeastOne® antifungal susceptibility testing of the new triazoles voriconazole, posaconazole and ravuconazole. *Journal of Clinical Microbiology* **42**: 4577-4580
- (12). Espinel-Ingroff, A., *et al.* (2004) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the NCCLS M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology* **42**: 718- 721
- (13) Linares,M, J., G. Charriel, F. Solís and M. Casal. (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology* **42**: 899-902
- (14) Espinel-Ingroff, A., M. Pfaller, S.A. Messer, C.C. Knapp, S. Killian, H.A. Norris, and M.A. Ghannoum. (1999) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology* **37**: 591-595
- (15) Barry, A.L., M.A. Pfaller, S.D. Brown, A. Espinel-Ingroff, M.A. Ghannoum, C.c. Knapp, R.P. Rennie, J.H. Rex, and M.G. Rinaldi. (2000) Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology* **38**: 3457-3459
- (16) Canton, E., *et al* (2005) .Sensititre YeastOne® Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of CLSI M27-A2 multicenter study . *Antimicrobial Agents and Chemotherapy*. **49**: 1604 -1607
- (17) Holiday, N.M... , *et al* (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne® susceptibility plate. *American Microbiology Abstract C-191*
- (18) Torres-Nabona, M., J. Guinea, J.Martinez_Alarcon, T. Pelaez, and E. Bouza. (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of NCCLS M38-A , Sensititre_YeastOne®. and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 - 1129
- (19) Espinel-Infgrroff, A. (2006) Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622

- (20) Patel, R., C. Mendrick, C.C. Knapp, R. Grist, and P.M. McNicholas. (2007) Clinical evaluation of the Sensititre YeastOne® plate for testing susceptibility of filamentous fungi to posaconazole. *Journal of Clinical Microbiology* **45**: 2000-2001
- (21) Alexander, B. D., T.C. Byrne, K.L. Smith, K.E. Hansen, K.J. Anstrom, J.R. Perfect and L. B. Reller. (2007) Comparative evaluation of Etest and Sensititre YeastOne® panels against the Clinical Laboratory Standards Institute M27-A2 reference broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology* **45**: 698-706
- (22). Patel, R., *et al.* (2007). Clinical evaluation of the Sensititre YeastOne plate for testing susceptibility of filamentous fungi to posiconazole. *Journal of Clinical Microbiology*. **45**: 2000-2001

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029-Yeast-ROW-IVD-I-V2.0
Date Updated: 06 August 2012

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LITERATUR

- (1) Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087.
- (2) Espinel-Ingroff, A., Pfaller, M., Messer, S.A., Knapp, C.C., Killian, S., Norris, H.A., Ghannoun M.A. (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology* **37**: 591-595.
- (3) Espinel-Ingroff, A., Knapp, C.C., Holiday, N., Killian, S. (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLSI M27-A microdilution reference method. *Clinical Microbiology and Infection* **8**: Suppl 1. Abstract 0125.
- (4) Davey, K.G., Szekely, A., Johnson, E.M., Warnock,D.W., (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Meletiadis, J., Moutin, J.W., Meis, F.G.M., Bouman, B.A., Verweij, P.E., and EuroFung Network. (2002). Comparison of the E Test and the Sensititre® Colorimetric Methods with the NCCLS Proposed Standard for Antifungal Susceptibility Testing of *Aspergillus* species. *Journal of Clinical Microbiology*. **40**: 2876-2885.
- (6) Espinel-Ingroff, A., etal. (1997). Multicenter Evaluation of Proposed Standardized Procedure for Antifungal Susceptibility Testing of Filamentous Fungi. *Journal of Clinical Microbiology* **35**: 139-143.

- (7) Pemán, J., et al. (2001). Comparison of the Sensititre YeastOne® Colorimetric method against reference M38-P for testing susceptibility of *Aspergillus* to Amphotericin B and Itraconazole. *Clinical Microbiology and Infection*, **7**: Suppl. 1. Abstract P694.
- (8) Linares, M., J. G. Charriel, F. Solís, F. Rodriguez, A. Ibarra and M. Casal. (2005) Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253
- (9) Kartsonis, N., et al. (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623
- (10) Pfaller, M, A., et al. (2004) Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases* **50** :113-117
- (11). Pfaller, M, A., A. Espinel-Ingroff and R.N. Jones. (2004) Clinical evaluation of the Sensititre YeastOne® antifungal susceptibility testing of the new triazoles voriconazole, posaconazole and ravuconazole. *Journal of Clinical Microbiology* **42**: 4577-4580
- (12). Espinel-Ingroff, A., et al. (2004) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the NCCLS M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology* **42**: 718- 721
- (13) Linares, M., J., G. Charriel, F. Solís and M. Casal. (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology* **42**: 899-902
- (14) Espinel-Ingroff, A., M. Pfaller, S.A. Messer, C.C. Knapp, S. Killian, H.A. Norris, and M.A. Ghannoum. (1999) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology* **37**: 591-595
- (15) Barry, A.L., M.A. Pfaller, S.D. Brown, A. Espinel-Ingroff, M.A. Ghannoum, C.c. Knapp, R.P. Rennie, J.H. Rex, and M.G. Rinaldi. (2000) Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology* **38**: 3457-3459
- (16) Canton, E., et al (2005) .Sensititre YeastOne® Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of NCCLS M27-A2 multicenter study . *Antimicrobial Agents and Chemotherapy*. **49**: 1604 -1607
- (17) Holiday, N.M., et al (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne® susceptibility plate. *American Microbiology Abstract C-191*
- (18) Torres-Nabona, M., J. Guinea, J.Martinez_Alarcon, T. Pelaez, and E. Bouza. (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of NCCLS M38-A ,

Sensititre_YeastOne[®] and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 - 1129

- (19) Espinel-Infgriff, A. (2006) Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622
- (20) Patel, R., C. Mendrick, C.C. Knapp, R. Grist, and P.M. McNicholas. (2007) Clinical evaluation of the Sensititre YeastOne[®] plate for testing susceptibility of filamentous fungi to posaconazole. *Journal of Clinical Microbiology* **45**: 2000-2001
- (21) Alexander, B. D., T.C. Byrne, K.L. Smith, K.E. Hansen, K.J. Anstrom, J.R. Perfect and L. B. Reller. (2007) Comparative evaluation of Etest and Sensititre YeastOne[®] panels against the Clinical Laboratory Standards Institute M27-A2 reference broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology* **45**: 698-706
- (22). Patel, R., etal. (2007). Clinical evaluation of the Sensititre YeastOne plate for testing susceptibility of filamentous fungi to posiconazole. *Journal of Clinical Microbiology*. **45**: 2000-2001

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029-Yeast-ROW-IVD-D-V2.0
Date Updated: 6 August 2012

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(1) Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087.

(2) Espinel-Ingroff, A., Pfaller, M., Messer, S.A., Knapp, C.C., Killian, S., Norris, H.A., Ghannoun M.A. (1999). Multicentre comparison of the Sensititre YeastOne[®] Colorimetric Antifungal Plate with the NCCLS, Institute M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology* **37**: 591-595.

(3) Espinel-Ingroff, A., Knapp, C.C., Holiday, N., Killian, S. (2002). Sensititre YeastOne[®] colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection* **8**: Suppl 1. Abstract 0125.

- (4) Davey, K.G., Szekely, A., Johnson, E.M., Warnock,D.W., (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Meletiadis, J., Moutin, J.W., Meis, F.G.M., Bouman, B.A., Verweij, P.E., and EuroFung Network. (2002). Comparison of the E Test and the Sensititre® Colorimetric Methods with the NCCLS Proposed Standard for Antifungal Susceptibility Testing of *Aspergillus* species. *Journal of Clinical Microbiology*. **40**: 2876-2885.
- (6) Espinel-Ingroff, A., *et al.* (1997). Multicenter Evaluation of Proposed Standardized Procedure for Antifungal Susceptibility Testing of Filamentous Fungi. *Journal of Clinical Microbiology* **35**: 139-143.
- (7) Pemán, J., *et al.* (2001). Comparison of the Sensititre YeastOne® Colorimetric method against reference M38-P for testing susceptibility of *Aspergillus* to Amphotericin B and Itraconazole. *Clinical Microbiology and Infection*, **7**: Suppl. 1. Abstract P694.
- (8) Linares,M, J., G. Charriel, F. Solís, F. Rodriguez, A. Ibarra and M. Casal. (2005) Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253
- (9) Kartsonis, N., *et al.* (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623
- (10) Pfaller, M, A., *et al.* (2004) Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases* **50** :113-117
- (11). Pfaller, M, A., A. Espinel-Ingroff and R.N. Jones. (2004) Clinical evaluation of the Sensititre YeastOne® antifungal susceptibility testing of the new triazoles voriconazole, posiconazole and raruconazole. *Journal of Clinical Microbiology* **42**: 4577-4580
- (12). Espinel-Ingroff, A., *et al.* (2004) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the NCCLS M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology* **42**: 718- 721
- (13) Linares,M, J., G. Charriel, F. Solís and M. Casal. (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology* **42**: 899-902
- (14) Espinel-Ingroff, A., M. Pfaller, S.A. Messer, C.C. Knapp, S. Killian, H.A. Norris, and M.A. Ghannoum. (1999) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology* **37**: 591-595
- (15) Barry, A.L., M.A. Pfaller, S.D. Brown, A. Espinel-Ingroff, M.A. Ghannoum, C.c. Knapp, R.P. Rennie, J.H. Rex, and M.G. Rinaldi. (2000) Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology* **38**: 3457-3459

- (16) Canton, E., *etal* (2005) .Sensititre YeastOne® Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of NCCLS M27-A2 multicenter study . *Antimicrobial Agents and Chemotherapy*. **49**: 1604 -1607
- (17) Holiday, N.M., *etal* (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne® susceptibility plate. *American Microbiology Abstract C-191*
- (18) Torres-Nabona, M., J. Guinea, J.Martinez_Alarcon, T. Pelaez, and E. Bouza. (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of NCCLS M38-A , Sensititre_YeastOne® and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 - 1129
- (19) Espinel-Infgriff, A. (2006) Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622
- (20) Patel, R., C. Mendrick, C.C. Knapp, R. Grist, and P.M. McNicholas. (2007) Clinical evaluation of the Sensititre YeastOne® plate for testing susceptibility of filamentous fungi to posaconazole. *Journal of Clinical Microbiology* **45**: 2000-2001
- (21) Alexander, B. D., T.C. Byrne, K.L. Smith, K.E. Hansen, K.J. Anstrom, J.R. Perfect and L. B. Reller. (2007) Comparative evaluation of Etest and Sensititre YeastOne® panels against the Clinical Laboratory Standards Institute M27-A2 refernce broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology* **45**: 698-706
- (22). Patel, R., *etal*. (2007). Clinical evaluation of the Sensititre YeastOne plate for testing susceptibility of filamentous fungi to posiconazole. *Journal of Clinical Microbiology*. **45**: 2000-2001

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029-Yeast-ROW-IVD-GR-V2.0
Date Updated: 06 August 2012

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BIBLIOGRAPHIE

- (1) Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087.
- (2) Espinel-Ingroff, A., Pfaller, M., Messer, S.A., Knapp, C.C., Killian, S., Norris, H.A., Ghannoun M.A. (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology* **37**: 591-595.
- (3) Espinel-Ingroff, A., Knapp, C.C., Holiday, N., Killian, S. (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection* **8**: Suppl 1. Abstract 0125.
- (4) Davey, K.G., Szekely, A., Johnson, E.M., Warnock,D.W., (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Meletiadis, J., Moutin, J.W., Meis, F.G.M., Bouman, B.A., Verweij, P.E., and EuroFung Network. (2002). Comparison of the E Test and the Sensititre® Colorimetric Methods with the NCCLS Proposed Standard for Antifungal Susceptibility Testing of *Aspergillus* species. *Journal of Clinical Microbiology*. **40**: 2876-2885.
- (6) Espinel-Ingroff, A., *et al.* (1997). Multicenter Evaluation of Proposed Standardized Procedure for Antifungal Susceptibility Testing of Filamentous Fungi. *Journal of Clinical Microbiology* **35**: 139-143.
- (7) Pemán, J., *et al.* (2001). Comparison of the Sensititre YeastOne® Colorimetric method against reference M38-P for testing susceptibility of *Aspergillus* to Amphotericin B and Itraconazole. *Clinical Microbiology and Infection*, **7**: Suppl. 1. Abstract P694.

- (8) Linares, M., J., G. Charriel, F. Solís, F. Rodriguez, A. Ibarra and M. Casal. (2005) Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253
- (9) Kartsonis, N., etal. (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623
- (10) Pfaller, M, A., etal. (2004) Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases* **50** :113-117
- (11). Pfaller, M, A., A. Espinel-Ingroff and R.N. Jones. (2004) Clinical evaluation of the Sensititre YeastOne® antifungal susceptibility testing of the new triazoles voriconazole, posaconazole and ravuconazole. *Journal of Clinical Microbiology* **42**: 4577-4580
- (12). Espinel-Ingroff, A., etal. (2004) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the NCCLS M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology* **42**: 718- 721
- (13) Linares,M, J., G. Charriel, F. Solís and M. Casal. (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology* **42**: 899-902
- (14) Espinel-Ingroff, A., M. Pfaller, S.A. Messer, C.C. Knapp, S. Killian, H.A. Norris, and M.A. Ghannoum. (1999) Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology* **37**: 591-595
- (15) Barry, A.L., M.A. Pfaller, S.D. Brown, A. Espinel-Ingroff, M.A. Ghannoum, C.c. Knapp, R.P. Rennie, J.H. Rex, and M.G. Rinaldi. (2000) Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology* **38**: 3457-3459
- (16) Canton, E., etal (2005) .Sensititre YeastOne® Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of NCCLS M27-A2 multicenter study . *Antimicrobial Agents and Chemotherapy*. **49**: 1604 -1607
- (17) Holiday, N.M..., etal (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne® susceptibility plate. *American Microbiology Abstract C-191*
- (18) Torres-Nabona, M., J. Guinea, J.Martinez_Alarcon, T. Pelaez, and E. Bouza. (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of NCCLS M38-A , Sensititre_YeastOne®. and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 - 1129
- (19) Espinel-Infroff, A. (2006) Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing

Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622

(20) Patel, R., C. Mendrick, C.C. Knapp, R. Grist, and P.M. McNicholas. (2007) Clinical evaluation of the Sensititre YeastOne® plate for testing susceptibility of filamentous fungi to posaconazole. *Journal of Clinical Microbiology* **45**: 2000-2001

(21) Alexander, B. D., T.C. Byrne, K.L. Smith, K.E. Hansen, K.J. Anstrom, J.R. Perfect and L. B. Reller. (2007) Comparative evaluation of Etest and Sensititre YeastOne® panels against the Clinical Laboratory Standards Institute M27-A2 reference broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology* **45**: 698-706

(22). Patel, R., *et al.* (2007). Clinical evaluation of the Sensititre YeastOne plate for testing susceptibility of filamentous fungi to posiconazole. *Journal of Clinical Microbiology*. **45**: 2000-2001

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029-Yeast-ROW-IVD-F-V2.0
Date Updated: 6 August 2012

Symbols/Símbolos/Simboli/Symbol/Symboles/σύμβολο/Simbolis/Symbole

LOT	GB Batch code DE Chargenbezeichnung ES Código de lote IT Codice del lotto FR Code du lot GR Κωδικός παρτίδας LI serijos numeri. CZ Šarže kadex PL Kod Partii RO lot număr PT Código do grupo	
REF	GB Catalogue number DE Bestellnummer ES Número de catálogo IT Numero di catalogo FR Référence du catalogue GR Αριθμός καταλόγου LI Reikia žinoti plokštelės kodą CZ Katalogovho číslovka PL Numer Katalogu RO Catalog număr PT Número de catálogo	
	GB Manufacturer DE Hersteller ES Fabricante IT Fabbricante FR Fabricant GR Κατασκευαστής LI Dirbt CZ Výrobce PL Producent RO Manufacturer PT Fabricante	
	GB Consult Instructions for Use- DE Gebrauchsanweisung beachten ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso FR Consulter les instructions d'utilisation GR Συμβουλευτείτε τις οδηγίες χρήσης LI Ieskant instrukciju CZ Projednat Instrukce Zauziti PL Przed użyciem zapoznaj się z instrukcją RO Consulta instrucțiuni pentru folos PT Consulte instruções para o uso	
	GB Temperature limitation DE Zulässiger Temperaturbereich ES Límite de temperatura IT Limiti di temperatura FR Limites de température GR Οπο θερμοκρασίας LI Temperatura reba CZ Teplota Limit PL Ograniczenie Temperatur RO Temperatură limitation PT Limitação de temperatura	
	GB Use By DE Verwendbar bis ES Fecha de caducidad IT Utilizzare entro FR Date de péremption GR Ημερομηνία λήξης LI Baigtis Data CZ Užití Za PL Użyć Do RO Folos By PT Uso perto	YYYYMMDD/YYYYYMM (MM=end of month) JJJJMMTT/JJJJmm (MM=Monatsende) aaaamddd/aaaamm (mm=fin del mes) AAAAMMGG/AAAAMM (MM=fine mese) AAAAMMJJ/AAAAMM (MM=fin de mois) EEEE-MM-HH / EEEE-MM (MM = τέλος του μήνα) KKKKMMDD RRRRMMDD RRRR/MM/DD AAAA/LL/ZZ YYYYMMDD/YYYYYMM (MM= fim do mês)
IVD	GB In Vitro Diagnostic Medical Device DE In Vitro Diagnostikum ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro FR Dispositif médical de diagnostic in vitro GR Ιν νιτρο διαγνωστικό ιατρικό βοήθημα LI In vitro diagnostikai CZ In vitro diagnostic PL Urządzenie Diagnostyczne in vitro RO Pentru diagnosticare in vitro PT In Vitro Dispositivo Médico Diagnóstico	

	GB Contains sufficient for <n> tests DE Ausreichend für <n> Ansätze ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per <n> saggi FR Contenu suffisant pour <n> tests GR Περιέχει επαρκή ποσατητά για <n> εξετάσεις LI Skaičius bandymas <n> CZ Obsahuje vhodný <n> testy PL Wystarczający dla <n> próbek RO contact sufficient pentru <n> tests PT Contem suficiente para <n> testes
	GB Plates with fluorescence substrates in wells DE Platten mit Fluoreszenzsubstraten in den Mulden ES Placa con substratos fluorescentes en los pocillos IT Pannelli con substrati fluorescenti nei pozzetti FR Plaquette avec substrats fluorescents dans les cuves GR Πλάκα με υποστρώματα φθορισμού σε υποδοχές LI substrato plokštelės CZ Fluorescenční substráty PL Płytki z substratami fluorescencyjnymi w studzienkach RO Fluorescentă Substratul PT Placas com as carcaças da fluorescência nos poços

SYMBOLS_V2.0